

BC-7500 CRP

Auto Hematology Analyzer

Principles

SF Cube method to count WBC, 6-part diff, NRBC,RET* and PLT-O*
(S: Scatter; F: Fluorescence; Cube: 3D analysis)
Focusing Flow-DC Impedance method for RBC and PLT
Cyanide free reagent for hemoglobin test
Latex Particle-Enhanced Immunonephelometry method for CRP test

Parameters

40 Reportable parameters(whole blood):

WBC, Lym%, Mon%, Neu%, Bas%, Eos%, IMG#, Lym#, Mon#, Neu#,
Eos#, Bas#, IMG#, RBC, HGB, HCT, MCV, MCH, MCHC, RDW-CV, RDW-SD,
NRBC#, NRBC%; PLT, MPV, PDW, PCT, P-LCR, P-LCC, RET%*, RET#*, RHE*,
IRF*, LFR*, MFR*, HFR*, IPF*, FR-CRP, CRP, hs-CRP

29 Research parameters(whole blood):

HFC#, HFC%, RBC-O*, PLT-O*, PLT-I, WBC-O*, WBC-D, TNC-D, IME%,
IME#, H-NR%, L-NR%, NLR, PLR, WBC-N, TNC-N, InR#, InR%o, Micro#,
Micro%, Macro#, Macro%, RPI*, H-IPF*, IPF#*, MRV*, FRC#*, FRC%*,
PDW-SD

7 Reportable parameters(body fluid):

WBC-BF, TC-BF#, MN#, MN%, PMN#, PMN%, RBC-BF

11 Research parameters (body fluid):

Eos-BF#, Eos-BF%, Neu-BF#,Neu-BF%, HF-BF#, HF-BF%, RBC-BF, LY-BF#,
LY-BF%, MO-BF#, MO-BF%

2 Histogram for RBC and PLT

3 Three-dimension scatter grams: DIFF, WNB, RET*

5 Two-dimension scatter grams: DIFF, WNB, RET*, RET-EXT*, PLT-O*

Modes

CBC, CBC+DIFF, CBC+DIFF+RET*, CBC+RET*, RET*

CBC+DIFF+CRP, CBC+DIFF+RET*+CRP, CRP

CBC+RET/PLT-8X*, CBC+DIFF+RET/PLT-8X*

Sample type	Sample volume
Whole Blood	≤ 85ul
Capillary Blood *	CRP mode: ≤ 20 μL CBC, CD mode: ≤ 27.5 μL CDR*, RET*, CD+CRP, and other modes ≤ 36.5 μL (+: When using dedicated capillary blood tubes provided by Mindray)
Predilute Sample	≤ 20ul
Body Fluid	≤ 20ul

Sample type	Analysis modes	Throughput
Whole Blood	CBC/CD	Up to 110 T/H
	CD+CRP	Up to 100 T/H
	CDR/CDR+CRP	Up to 65 T/H
Capillary Blood *	CBC/CD/CD+CRP	Up to 70 T/H
	CDR/CDR+CRP	Up to 45 T/H
(*: When using dedicated capillary blood tubes provided by Mindray)		

Performance

Parameter	Linearity Range	Precision	Carryover
WBC	0~500.00×10 ⁹ /L	≤ 2.5% (≥ 4.51×10 ⁹ /L)	≤1.0%
RBC	0.00~8.60×10 ¹² /L	≤ 1.5% (≥3.50×10 ¹² /L)	≤1.0%
HGB	0~260 g/L	≤ 1.0% (110 ~ 180) /L	≤1.0%
HCT	0.0~75.0%	≤ 1.5% (30.0 ~ 50.0)%	≤1.0%
PLT	0~5000×10 ⁹ /L	≤ 4.0% (≥100×10 ⁹ /L)	≤1.0%
RET#	(0.0000~0.8000)×10 ¹² /L	≤ 15%(RBC≥ 3.00×10 ¹² /L RET%:1.00%~4.00%)	/
FR-CRP	0.20mg/L~320.00mg/L	≤0.50mg/L(SD) (≤10mg/L) ≤4% (CV) (>10mg/L)	≤1.0%

Loading capacity

Up to 50 sample tubes

Data storage capacity

Up to 150,000 results including numeric and graphical information

Operating environment

Ambient temperature: 15℃ ~ 32℃

Relative humidity: 30% ~ 85%

Atmospheric pressure: 70.0 kPa ~ 106.0 kPa

*Only available in Reticulocyte configuration

BC-7500 CRP

Auto Hematology Analyzer

Empowering Pediatric Care




Capillary blood autoloading and mixing



Capillary blood rack



Dedicated capillary blood tube



Biotic mixing component

- Eccentric axis and continuously start and stop rotation
- Alternation between clockwise and counterclockwise

- Fully exchange of upper and lower blood samples ensures the best mixing performance of capillary blood

High speed CBC+CRP integrated analyzer

100T/H
CD+CRP

Aspirating needle, Light source, Sample pool, LC hemolytic reagent, Signal receiver, Whole process temperature control, Latex

Constant temperature **FAST**

- + Preheat the reagent to reach the best reaction temperature faster

Sufficient mixing **FAST**

- + Control the speed, angle and time when adding reagent to enhance hemolysis and immune reaction

Efficient lyse reagent **FAST**

- + Efficient and specific multi-functional lyse reagent targets at the rigid components of cell membrane skeleton for fast dissolving

Whole blood

Plasma, PLT, WBC, RBC

Serum

CRP antigens, Serum matrix

$$CRP_{Serum} = \frac{CRP_{Whole\ blood}}{1 - BCT\%}$$



BCT: hematocrit of RBC, WBC and PLT

Stable reaction **FAST**

- + Pre-reaction hemolytic reagent
- + High signal latex reagent
- + Fixed site reaction
- + Early rapid reaction segment detection


Reliable result **ACCURATE**

- + The function of CRP modification ensures that the result of CRP_{Serum} and $CRP_{whole\ blood}$ are consistent.




DIFF Channel

- WBC 6-part differential results including immature granulocyte.
- Provides research parameters such as HFC (Blast & Atypical Lymphocyte), InR (information about malaria) and flags for Band, NRBC, PLT clump and Atypical Lymphocyte.
- 3-times counting of low WBC ensures more reliable counting and differentiation result




WNB Channel

- Provides WBC, NRBC and Basophils results.
- NRBC in every CBC test helps to obtain reliable WBC result for neonatal sample
- It also helps the diagnosis for hemolytic anemia, hematopoietic diseases and reduce the review rate.



RET Channel

- Provides RET related parameters (i.e. RET#/RET%/LFR/MFR/HFR/RHE) which helps anemia diagnosis
- Fluorescent PLT result helps to remove the interference of fragmented RBC/PLT clump
- /large PLT, etc.
- Automatic PLT-O* 8 times extended counting helps to obtain reliable PLT result even in low count.



BF detection

- Can also detect 4 types of body fluids.